



**LABORATORY ALLIANCE**



of Central New York, LLC

laboratoryalliance.com


## Labeling Placement Change Announced

Please note the following changes to label placement on these BacT/ALERT® culture bottles: aerobic SA and anaerobic SN culture bottles.

Place patient labels in the area indicated in **green ONLY** for ideal Blood Level Detection performance. **Red areas are not appropriate for label placement.**



The label still includes all of the current elements, a fill line and a 2D barcode.



2D Barcode

Fill-to mark  
(white zone  
etched)

For more information, email Microbiology Manager Russell Rawling at [russellrawling@lacny.com](mailto:russellrawling@lacny.com) or call our Customer Service Department at 315-461-3008.