



LABORATORY ALLIANCE


of Central New York, LLC

laboratoryalliance.com

Labeling Placement Change Announced

Please note the following changes to label placement on this BacT/ALERT® culture bottle: pediatric PF Plus culture bottle.

Place patient labels in the area indicated in **green** ONLY for ideal Blood Level Detection performance. **Red areas are not appropriate for label placement.**



2D Barcode

The label still includes all of the current elements and a 2D barcode.



For more information, email Microbiology Manager Russell Rawling at russellrawling@lacny.com or call our Customer Service Department at 315-461-3008.

Laboratory Alliance of CNY 1/2020