

## Smart Transport Media (STM), made by MedSchenker, can be used in place of Universal Transport Media (UTM) for the following tests:

- COVID-19 by PCR performed at Lab Alliance, University Hospital, ARUP or LabCorp
- Flu A, Flu B and RSV by PCR
- Herpes Simplex Virus 1 & 2/Varicella Zoster Virus, Molecular
- Bordetella pertussis by PCR
- Respiratory panel: SJH and lab-order-only Respiratory Panel to University Hospital



NOTE: STM can be stored at Room Temperature (20-25°C) or Refrigerated (2-8°C).

Please contact the Microbiology Department if you have any questions. (315) 410-7067