



LABORATORY ALLIANCE
of Central New York, LLC

Viral Transport Media Substitute for Universal Transport Media (UTM)

Note: Must be stored at 2 – 8⁰C

VTM can be used for these tests requiring UTM:

- COVID-19 by PCR
- Flu A, Flu B and RSV by PCR
- Herpes Simplex Virus 1 & 2/Varicella Zoster Virus Molecular
- Bordetella pertussis PCR
- Respiratory panel: SJH and Lab Order Only
Respiratory Panel to SUNY

